




176-LED-GM

GM Hammered Grey


176-LED-BM

BM Hammered White



176-LED

ACCESSORIES



S11

S12MED

\$15

RL

BASEPE

Power	15.5 W
Supply	110–240 V
Illuminance at 50cm	1.770 lux
Lighting Intensity Adjustment	0,5 – 100 %
Colour Temperature	5.300 K
Reflector's length	40 cm
Diameter of the light field	700 x 650 mm
Product's life	50.000 hours
Structure	Reflector + pantograph arm with spring
IP Protection	IP20
IK Protection	IK05
Class	III
Cabling	Plug-top transformer with AWG 2X21; interchangeable plugs
Fixing	S11 / S12MED / S15 clamp / RL / BASEPE
Colours	Hammered Grey, Hammered White
Order Code	176-LED-GM, 176-LED-BM

